

Mode	Mode in binary	Activity	Dual Activity	3rd Axis	Steps	Heart Rate	Lux	Incline
0	000000	X						
1	000001	X			X			
2	000010	X				X		
3	000011	X			X	X		
4	000100	X	X					
5	000101	X	X		X			
6	000110	X	X			X		
7	000111	X	X		X	X		
8	001000	X		X				
9	001001	X		X	X			
10	001010	X		X		X		
11	001011	X		X	X	X		
12	001100	X	X	X				
13	001101	X	X	X	X			
14	001110	X	X	X		X		
15	001111	X	X	X	X	X		
16	010000	X					X	
17	010001	X			X		X	
18	010010	X				X	X	
19	010011	X			X	X	X	
20	010100	X	X				X	
21	010101	X	X		X		X	
22	010110	X	X			X	X	
23	010111	X	X		X	X	X	
24	011000	X		X			X	
25	011001	X		X	X		X	
26	011010	X		X		X	X	
27	011011	X		X	X	X	X	
28	011100	X	X	X			X	
29	011101	X	X	X	X		X	
30	011110	X	X	X		X	X	
31	011111	X	X	X	X	X	X	
32	100000	X						X
33	100001	X			X			X
34	100010	X				X		X
35	100011	X			X	X		X
36	100100	X	X					X
37	100101	X	X		X			X
38	100110	X	X			X		X
39	100111	X	X		X	X		X
40	101000	X		X				X
41	101001	X		X	X			X
42	101010	X		X		X		X
43	101011	X		X	X	X		X
44	101100	X	X	X				X
45	101101	X	X	X	X			X
46	101110	X	X	X		X		X
47	101111	X	X	X	X	X		X
48	110000	X					X	X
49	110001	X			X		X	X
50	110010	X				X	X	X
51	110011	X			X	X	X	X
52	110100	X	X				X	X
53	110101	X	X		X		X	X
54	110110	X	X			X	X	X
55	110111	X	X		X	X	X	X
56	111000	X		X			X	X
57	111001	X		X	X		X	X
58	111010	X		X		X	X	X
59	111011	X		X	X	X	X	X
60	111100	X	X	X			X	X
61	111101	X	X	X	X		X	X
62	111110	X	X	X		X	X	X
63	111111	X	X	X	X	X	X	X